



Enteral nutrition: Administration of nutritional formula - Button - By pump

- Anticipate
- Prepare
- Proceed
- Prepare the material
- Position your child
- Check the location of the button
- Irrigate the button extension set
- Connect the extension set to the button
- Irrigate the button
- Pour the formula into the bag
- Hang the bag
- Prepare the pump
- Empty the air out of the bag tubing
- Connect the tubing to the button extension set
- Start the pump
- Administer the formula
- Stop the pump
- Disconnect the tubing from the button extension set
- Irrigate the button
- Disconnect the button extension set
- Keep your child in the same position after the session
- Clean and store equipment
- Recover