

Intravenous antibiotic therapy: Administration of the intravenous antibiotics by ambulatory infusion pump: Change of bag AND tubing

pump: Change of bag AND tubing
Anticipate
Prepare
Proceed
Prepare the drug
Prepare the materials
Position your child
Stop the pump
Remove the cassette from the pump
Install a new bag and cassette tubing
Remove the air from the new tubing
Disconnect the old tubing from the needleless connector
Monitor the catheter site
Irrigate the catheter with saline solution, if necessary
Connect the new tubing to the needleless connector of the catheter
Restart the pump
Secure the pump and medication
Recover