



Intravenous antibiotic therapy: Administration of the intravenous antibiotics by ambulatory infusion pump – Change of bag AND tubing
SUMMARY OF STEPS

Steps	Notes
Anticipate	
Prepare	
Proceed	
Prepare the drug	
Prepare the materials	
Position your child	
Stop the pump	
Remove the cassette from the pump	
Install a new bag and cassette tubing	
Remove the air from the new tubing	
Disconnect the old tubing from the needleless connector	



Intravenous antibiotic therapy: Administration of the intravenous antibiotics by elastomeric infusion pump
Change of bag AND tubing

SUMMARY OF STEPS

Steps	Notes
Monitor the catheter site	
Irrigate the catheter with saline solution, if necessary	
Connect the new tubing to the needleless connector of the catheter	
Restart the pump	
Secure the pump and medication	
Recover	