



Intravenous antibiotic therapy: Administration of the intravenous antibiotics by ambulatory infusion pump – Change of cassette reservoir AND tubing
SUMMARY OF STEPS

| Steps | Notes |
|---|-------|
| Anticipate | |
| Prepare | |
| Proceed | |
| Prepare the drug | |
| Prepare the materials | |
| Position your child | |
| Stop the pump | |
| Remove the empty cassette reservoir from the pump | |
| Install a new cassette reservoir and new extension tubing | |
| Remove the air from the new tubing | |
| Disconnect the old tubing from the needleless connector | |
| Monitor the catheter site | |

| Steps | Notes |
|--|-------|
| Irrigate the catheter with saline solution, if necessary | |
| Connect the new tubing to the needleless connector of the catheter | |
| Restart the pump | |
| Secure the pump and medication | |
| Recover | |