



**Intravenous antibiotic therapy: Administration of the intravenous antibiotics by ambulatory infusion pump – Change of cassette reservoir AND tubing**  
**SUMMARY OF STEPS**

Steps	Notes
Anticipate	
Prepare	
Proceed	
Prepare the drug	
Prepare the materials	
Position your child	
Stop the pump	
Remove the empty cassette reservoir from the pump	
Install a new cassette reservoir and new extension tubing	
Remove the air from the new tubing	
Disconnect the old tubing from the needleless connector	
Monitor the catheter site	
Irrigate the catheter with saline solution, if necessary	
Connect the new tubing to the needleless connector of the catheter	
Restart the pump	
Secure the pump and medication	
Recover	