



## Peritoneal dialysis: Assembly of the cyclor

- Anticipate
- Prepare
- Proceed
- Turn on the cyclor
- Check the dialysate bags
- Add medication to the bags, if necessary
- Prepare the appropriate dialysate bags
  - Dianeal
  - Extraneal
  - Physioneal 2.5liter: break the cannula
  - Physioneal 5 liter: open the long and the short seal and check the contour of the bag
- Place the bag on the heater
- Prepare the cyclor set with cassette
- Prepare drainage bag and specimen bag (if needed)
- Start the cyclor
- Open the cyclor door
- Install the cassette of the cyclor set
- Close the cyclor door
- Install the organizer of the cyclor set
- Disinfect your hands



- Connect the specimen bag (if needed)
- Connect the drainage bag
- Prepare the connection of the bags
- Disinfect your hands
- Connect the bag to the heater (taking account of bag specifics)
- Connect the dialysate bags
- Connect the last injection bag (if needed)
- Check the connections
- Place the bags
- Priming the lines
- Check the dialysate level in the patient line
- Connect your child to the cyclor, if desired
- Recover