



Enterostomy care: Changing the collection device: 1 piece

- Anticipate
- Prepare
- Proceed
- Position your child
- Empty the bag
- Remove the collection device
- Clean the stoma and skin around it
- Examine the stoma and the condition of the skin around it
- Place a temporary absorbent protection over the stoma
- Cut the opening in the skin barrier, if necessary
- Even the cut opening of the skin barrier
- Warm the skin barrier
- Select the placement of the bag
- Remove the plastic covering or the protective paper from the skin barrier
- Protect the skin around the stoma
- Place the skin barrier around the stoma
- Empty air out of the bag
- Close the bag
- Check that the skin barrier is completely sealed to the skin
- Recover