

Ventilatory assistance: Invasive ventilation

- Anticipate
- Prepare
- Proceed
- Prepare the device
- Assemble the ventilation circuit
- Prepare the heated humidifier, if being used
- Prepare the artificial nose, if being used
- Position and prepare your child
- Inflate the cuffed tracheal cannula, if present
- Turn on the device
- Connect the circuit to the tracheal cannula
- Continue and supervise ventilator assistance
- Discontinue ventilation and disconnect the circuit from the tracheal cannula
- Turn off the device
- Deflate the cuffed tracheal cannula, if present
- Replace the speaking valve or the artificial nose, if needed
- Clean and store the equipment
- Recover