



Ventilatory assistance: Manual ventilation by tracheostomy

SUMMARY OF STEPS

Steps	Notes
Anticipate	
Prepare	
Proceed	
Prepare and verify the ventilator bag	
Position your child	
Prepare your child	
Inflate the cuffed tracheal cannula, if present	
Connect the ventilatory bag to the tracheal cannula	
Ventilate with the ventilator bag	
Disconnect the ventilator bag from the tracheal cannula	
Deflate the cuffed tracheal cannula, if present	
Replace the speaking valve or the artificial nose, if needed	
Restart invasive ventilation, if needed	
Clean and store the equipment	
Recover	